

# Praktikum Visual Computing

## Object Detection for Chest X-Ray Imaging

By Yihao Li  
Tutor: Sanner Autoine

# Chest ImaGenome Dataset

- a subset of the MIMIC-CXR dataset
- "golden dataset": manually annotated
  - 256 combinations of relation annotations

patient_id		Region	sentence
15184836	00046130-fd952ef0-57f2948d-491a16b4-5db3a18c.dcm	['cardiac silhouette', 'left hilar structures', 'mediastinum', 'right hilar structures', 'upper mediastinum']	The cardiomediastinal and hilar contours are within normal limits.
15184836	00046130-fd952ef0-57f2948d-491a16b4-5db3a18c.dcm	['left lung', 'right lung']	The lung fields are clear.
15184836	00046130-fd952ef0-57f2948d-491a16b4-5db3a18c.dcm	['left apical zone', 'left lung', 'right apical zone', 'right lung']	There is no pneumothorax, fracture or dislocation.
12930467	005043e2-a4e25d1d-aae26631-732a2db0-38412248.dcm	['left hilar structures', 'left lower lung zone', 'left lung', 'right hilar structures', 'right lower lung zone', 'right lung']	lower lungs are grossly clear, though there is vascular crowding


- over 670, 000 localized comparison relations ( improved, worsened, or no change) between the anatomical locations

# Chest ImaGenome Dataset

---

Average object size: 648 x 586


Average image size: 2675 x 2794




# Chest ImaGenome Dataset

## Object Size Biggest 3 :


Left lung: 972 x 1809



Right lung: 958 x 1794




Cardiac silhouette: 1070 x 777




# Chest ImaGenome Dataset

## Object Size Smallest 3 :


Cavoatrial junction: 314 x 224



Aortic arch: 252 x 226



Carina: 121 x 122



# Detection Method

Stages	Detection Method	Advantages	Disadvantages
One-stage	YOLO	<ul style="list-style-type: none"><li>• Transform object detection task into a regression problem</li><li>• Highly speed up the detection</li></ul>	<ul style="list-style-type: none"><li>• Position is not accurate</li><li>• Detection effect of small and dense instances is not efficient</li></ul>
	SSD	<ul style="list-style-type: none"><li>• Use feature layers of different scales extracted from the backbone</li><li>• Speed of detection is fast</li></ul>	<ul style="list-style-type: none"><li>• Features for smaller targets may be lost</li></ul>
Two-stage	R-CNN	<ul style="list-style-type: none"><li>• Extract and learn features from CNNs automatically</li><li>• Compute the neural network features on each of regions of interest</li></ul>	<ul style="list-style-type: none"><li>• Long time to acquire regional targets</li><li>• Feature extraction is complex</li></ul>
	Fast R-CNN	<ul style="list-style-type: none"><li>• Run the neural network once on the whole image</li><li>• Training speed is significantly enhanced</li></ul>	<ul style="list-style-type: none"><li>• Restricted by selective search algorithms</li></ul>
	Faster R-CNN	<ul style="list-style-type: none"><li>• Region Proposal Network(RPN) instead of selective search algorithms</li><li>• Improve detection speed</li></ul>	<ul style="list-style-type: none"><li>• Not satisfy the requirement of real-time detection due to two stage structure</li></ul>

## Mmdetection:

- Maintained
- Framework base: Pytorch
- Pre-trained models:
  - Fast R-CNN
  - Faster R-CNN
  - SSD
  - RetinaNet
  - YOLOv3
  - other state-of-the-art detection models


## COCO format

```
{
  "annotations": [
    {
      "id": 120,
      "image_id": 4,
      "category_id": 4,
      "area": 74529.0,
      "bbox": [245.0, 1759.0, 273.0, 273.0],
      "iscrowd": 0
    }
  ],
  "categories": [
    {
      "id": 1,
      "name": "left lung"
    }
  ],
  "images": [
    {
      "id": 4,
      "width": 3056,
      "height": 2544,
      "file_name": "image1.jpg"
    }
  ]
}
```

- Split data:8:1:1
- Data augmentation:
  - translation: shift small distance
  - rotation: less than 10 degrees




# Evaluations metrics



$$\text{Precision} = \frac{TP}{TP + FP} = \frac{TP}{\text{all detections}}$$



mAP: average of the precision values at different IoU thresholds

# Result



# Result

		Precision @[IoU = 0.50]	Precision @[IoU = 0.75]	Precision @[IoU = 0.50 : 0.95]
Faster_rcnn	Translation	88,8%	54,1%	52,4%
	Rotation	87,6%	56,5%	52,8%
	Translation + Rotation	87,9%	57,8%	53,9%
Yolov3	Translation	84,4%	40,5%	43,6%
	Rotation	84,2%	34,9%	40,7%
	Translation + Rotation	84,4%	39,4%	42,9%
SSD	Translation	87%	50,4%	49,2%
	Rotation	87,1%	47,2%	47,4%
	Translation + Rotation	87,1%	49,7%	48,3%





# Left & Right organs

	Faster_rcnn		Yolo		SSD	
	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]
left hemidiaphragm	76%	43,1%	63,1%	25,7%	72,8%	35,3%
right hemidiaphragm	89,5%	57,2%	86,4%	39,3%	85,6%	46,9%
left lung	100%	78,1%	100%	62,8%	100%	61,3%
right lung	99%	80,0%	95,7%	62%	99,8%	65,6%
left hilar structures	94,3%	54,9%	92,2%	51%	97,3%	53,9%
right hilar structures	92,5%	56%	92,5%	49,9%	96%	52,8%
left apical zone	96,7%	64,9%	97,6%	45,1%	97,6%	61,5%
right apical zone	95,6%	58,8%	89%	43,6%	92,5%	56,3%
left clavicle	85%	53,5%	76,9%	29,4%	85,7%	46%
right clavicle	86,9%	51,1%	75,1%	26,4%	85,1%	46,9%
left costophrenic angle	62,1%	28,8%	66,4%	32,3%	65,5%	29,3%
right costophrenic angle	78,6%	37,4%	82,3%	40,3%	81,2%	39,5%

# Big & Small organs

	Faster_rcnn		Yolo		SSD	
	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]	Precision @[IoU = 0.50]	Precision @[IoU = 0.50 : 0.95]
left lung	100%	78,1%	100%	62,8%	100%	61,3%
right lung	99%	80,0%	95,7%	62%	99,8%	65,6%
cardiac silhouette	98,9%	63%	96,1%	53,3%	96,8%	54,4%
left costophrenic angle	62,1%	28,8%	66,4%	32,3%	65,5%	29,3%
carina	56,1%	16,8%	44,1%	11,9%	38,8%	11,2%
cavoatrial junction	53,7%	16,1%	57,2%	17,2%	59,8%	18,4%

# Predictions vs. Ground Truth



Left lung

Right lung

Cardiac silhouette


Cavoatrial junction

Aortic arch

Carina






# Left & Right Lung





# Cavoatrial junction









Ground Truth

Faster RCNN

YOLO


SSD







# Thanks!

learning\_rate




acc  
tag: train/acc






acc\_val  
tag: val/acc\_val



loss\_val  
tag: val/loss\_val





```

Average Precision (AP) @[ IoU=0.50:0.50 | area= all | maxDets=100 ] = 0.875
Average Precision (AP) @[ IoU=0.50 | area= all | maxDets=1000 ] = -1.000
Average Precision (AP) @[ IoU=0.75 | area= all | maxDets=1000 ] = -1.000
Average Precision (AP) @[ IoU=0.50:0.50 | area= small | maxDets=1000 ] = 0.000
Average Precision (AP) @[ IoU=0.50:0.50 | area=medium | maxDets=1000 ] = 0.078
Average Precision (AP) @[ IoU=0.50:0.50 | area= large | maxDets=1000 ] = 0.878
Average Recall (AR) @[ IoU=0.50:0.50 | area= all | maxDets=100 ] = 0.900
Average Recall (AR) @[ IoU=0.50:0.50 | area= all | maxDets=300 ] = 0.900
Average Recall (AR) @[ IoU=0.50:0.50 | area= all | maxDets=1000 ] = 0.900
Average Recall (AR) @[ IoU=0.50:0.50 | area= small | maxDets=1000 ] = 0.000
Average Recall (AR) @[ IoU=0.50:0.50 | area=medium | maxDets=1000 ] = 0.090
Average Recall (AR) @[ IoU=0.50:0.50 | area= large | maxDets=1000 ] = 0.904

```

category	AP	category	AP	category	AP
right lung	1.000	left lung	1.000	left costophrenic angle	0.656
trachea	0.927	right costophrenic angle	0.776	cavoatrial junction	0.524
left mid lung zone	0.900	left clavicle	0.856	right hemidiaphragm	0.881
aortic arch	0.895	left hemidiaphragm	0.799	right apical zone	0.940
carina	0.470	left hilar structures	0.966	right hilar structures	0.933
right clavicle	0.909	svc	0.870	upper mediastinum	0.973
right atrium	0.834	right mid lung zone	0.907	left apical zone	0.968
left upper lung zone	0.999	right upper lung zone	0.969	left lower lung zone	0.897
right lower lung zone	0.917	cardiac silhouette	0.977	None	None

```

Average Precision (AP) @[ IoU=0.50:0.95 | area= all | maxDets=100 ] = 0.550
Average Precision (AP) @[ IoU=0.50 | area= all | maxDets=1000 ] = 0.875
Average Precision (AP) @[ IoU=0.75 | area= all | maxDets=1000 ] = 0.586
Average Precision (AP) @[ IoU=0.50:0.95 | area= small | maxDets=1000 ] = 0.000
Average Precision (AP) @[ IoU=0.50:0.95 | area=medium | maxDets=1000 ] = 0.023
Average Precision (AP) @[ IoU=0.50:0.95 | area= large | maxDets=1000 ] = 0.551
Average Recall (AR) @[ IoU=0.50:0.95 | area= all | maxDets=100 ] = 0.612
Average Recall (AR) @[ IoU=0.50:0.95 | area= all | maxDets=300 ] = 0.612
Average Recall (AR) @[ IoU=0.50:0.95 | area= all | maxDets=1000 ] = 0.612
Average Recall (AR) @[ IoU=0.50:0.95 | area= small | maxDets=1000 ] = 0.000
Average Recall (AR) @[ IoU=0.50:0.95 | area=medium | maxDets=1000 ] = 0.031
Average Recall (AR) @[ IoU=0.50:0.95 | area= large | maxDets=1000 ] = 0.614

```

category	AP	category	AP	category	AP
right lung	0.822	left lung	0.800	left costophrenic angle	0.324
trachea	0.506	right costophrenic angle	0.391	cavoatrial junction	0.172
left mid lung zone	0.584	left clavicle	0.549	right hemidiaphragm	0.585
aortic arch	0.515	left hemidiaphragm	0.435	right apical zone	0.631
carina	0.156	left hilar structures	0.563	right hilar structures	0.584
right clavicle	0.549	svc	0.472	upper mediastinum	0.601
right atrium	0.381	right mid lung zone	0.569	left apical zone	0.658
left upper lung zone	0.787	right upper lung zone	0.748	left lower lung zone	0.615
right lower lung zone	0.653	cardiac silhouette	0.642	None	None